



# Cell Strainer 細胞過篩器



## pluriStrainer®

一般過篩器 (適合一般樣品及多種應用)

- pluriStrainer®共有 **19** 種孔徑選擇：  
1, 5, 10, 15, 20, 30, 40, 50, 60, 70, 85, 100, 150, 200, 300, 400, 500, 750及1000µm · 滿足你的過篩需求！
- PET (Polyethylenterephthalat)-mesh (過篩膜材質)  
(其中40/ 70/ 100µm 三種孔徑另有nylon-mesh選擇)
- 可直接放於50ml離心管上過篩使用

德國  
製造

多種孔徑  
選擇

應用  
多元

不管過篩什麼 都超好用！

- 從檢體得到單一細胞
- 消化組織得到單一細胞或細胞分離
- 細胞,植物種子,寄生蟲卵等樣品過濾或分離
- Flow上機前樣品準備
- 微珠細胞分選
- 用於過篩黏稠的液體樣品  
(結合使用連接器環和注射器)

### 勝 pluriStrainer 多樣化的應用

其他  
品牌

#### 正放過篩

放到50ml離心管上·  
即可過濾樣品



#### 過濾完反轉

後·放到50ml離心  
管上·即可沖洗下  
較大的樣品



堆疊不同孔徑過濾  
器·即可同時過濾  
不同大小樣品



#### 加上漏斗

放到50ml離心管上·  
即可增加過濾樣品  
體積·最多達24ml



加上連接環·藉  
由開關可控制流速·  
適合用於研磨組織  
取細胞



加上連接環與  
注射器·  
可造成強制性的低壓  
來過濾較黏稠樣品



## pluriStrainer® Mini

迷你小體積過篩器 (適合小量體積的樣品)

- 共 **7** 種孔徑：5, 10, 20, 40, 70 100 and 200µm
- 適用放於 1.5 /2.0 ml 反應管, 15 ml 離心管, FACS™ tubes, 冷凍小管, 24 及 48 well plate等容器



pluriStrainer® Mini  
都適用這些容器過篩

品項 / 規格	貨號
<b>pluriStraine (PET-Mesh) 一般過篩器</b>	
pluriStrainer Set 1, 32 pcs, sterile, 小孔徑組合包 4 x (1µm, 5µm, 10µm, 15µm, 20µm and 30µm) pluriStrainer and 8 Connector Rings	43-50001-99
pluriStrainer Set 2, 24 pcs, sterile, 中孔徑組合包 4 x (40µm, 60µm, 70µm, 85µm, 100µm and 200µm).	43-50000-98
pluriStrainer Set 3, 24 pcs, sterile, 大孔徑組合包 4 x (70µm, 100µm, 200µm, 300µm, 400µm and 500µm).	43-50000-97
pluriStrainer 1µm, 25 pcs, sterile + 25pcs Connector Ring, sterile	43-50001-03
pluriStrainer 1µm, 25 pcs, sterile	43-50001-13
pluriStrainer 1µm, 50 pcs - non-sterile	43-50001-50
pluriStrainer 5µm, 25 pcs, sterile + 25pcs Connector Ring, sterile	43-50005-03
pluriStrainer 5µm, 25 pcs, sterile	43-50005-13
pluriStrainer 5µm, 50 pcs - non-sterile	43-50005-50
pluriStrainer 10µm, 25 pcs, sterile	43-50010-03
pluriStrainer 10µm, 50 pcs - non-sterile	43-50010-50
pluriStrainer 15µm, 25 pcs, sterile	43-50015-03
pluriStrainer 15µm, 50 pcs - non-sterile	43-50015-50
pluriStrainer 20µm, 25 pcs, sterile	43-50020-03
pluriStrainer 20µm, 50 pcs - non-sterile	43-50020-50
pluriStrainer 30µm, 25 pcs, sterile	43-50030-03
pluriStrainer 30µm, 50 pcs - non-sterile	43-50030-50
pluriStrainer 40µm, 25 pcs, sterile	43-50040-03
pluriStrainer 40µm, 50 pcs - non-sterile	43-50040-50
pluriStrainer 40µm, 50 pcs, sterile	43-50040-51
PluriStrainer® 50 µm, 25 pcs. – sterile	43-50050-03
PluriStrainer® 50 µm, 50 pcs. - non-sterile	43-50050-50
pluriStrainer 60µm, 25 pcs, sterile	43-50060-03
pluriStrainer 60µm, 50 pcs - non-sterile	43-50060-50
pluriStrainer 70µm, 25 pcs, sterile	43-50070-03
pluriStrainer 70µm, 50 pcs - non-sterile	43-50070-50
pluriStrainer 70µm, 50 pcs, sterile	43-50070-51
pluriStrainer 85µm, 25 pcs, sterile	43-50085-03
pluriStrainer 85µm, 50 pcs - non-sterile	43-50085-50
pluriStrainer 100µm, 25 pcs, sterile	43-50100-03
pluriStrainer 100µm, 50 pcs - non-sterile	43-50100-50
pluriStrainer 100µm, 50 pcs, sterile	43-50100-51
pluriStrainer 150 µm, 25 pcs. – sterile	43-50150-03
pluriStrainer 150 µm, 50 pcs. - non-sterile	43-50150-50
pluriStrainer 200µm, 25 pcs, sterile	43-50200-03
pluriStrainer 200µm, 50 pcs - non-sterile	43-50200-50
pluriStrainer 300µm, 25 pcs, sterile	43-50300-03
pluriStrainer 300µm, 50 pcs - non-sterile	43-50300-50
pluriStrainer 400µm, 25 pcs, sterile	43-50400-03
pluriStrainer 400µm, 50 pcs - non-sterile	43-50400-50
pluriStrainer 500µm, 25 pcs, sterile	43-50500-03
pluriStrainer 500µm, 50 pcs - non-sterile	43-50500-50
pluriStrainer 750 µm, 25 pcs. – sterile	43-50750-03
pluriStrainer 750 µm, 50 pcs. - non-sterile	43-50750-50
pluriStrainer 1000 µm, 25 pcs. – sterile	43-51000-03
pluriStrainer 1000 µm, 50 pcs. - non-sterile	43-51000-50

品項 / 規格	貨號
<b>pluriStraine (Nylon-Mesh) 一般過篩器</b>	
pluriStrainer 40µm, Nylon-Mesh, 25 pcs, sterile	43-57040-03
pluriStrainer 40µm, Nylon-Mesh, 50 pcs, non-sterile	43-57040-50
pluriStrainer 40µm, Nylon-Mesh, 50 pcs, sterile	43-57040-51
pluriStrainer 70µm, Nylon-Mesh, 25 pcs, sterile	43-57070-03
pluriStrainer 70µm, Nylon-Mesh, 50 pcs, non-sterile	43-57070-50
pluriStrainer 70µm, Nylon-Mesh, 50 pcs, sterile	43-57070-51
pluriStrainer 100µm, Nylon-Mesh, 25 pcs, sterile	43-57100-03
pluriStrainer 100µm, Nylon-Mesh, 50 pcs, non-sterile	43-57100-50
pluriStrainer 100µm, Nylon-Mesh, 50 pcs, sterile	43-57100-51
<b>pluriStrainer Mini (PET-Mesh) 迷你小體積過篩器</b>	
pluriStrainer Mini, 5µm, 25pcs sterile in bag	43-10005-40
pluriStrainer Mini, 5µm, 50pcs, 25pcs sterile in bag	43-10005-50
pluriStrainer Mini, 5µm, 100pcs, 25pcs sterile in bag	43-10005-60
pluriStrainer Mini, 5µm, 500pcs, 25pcs sterile in bag	43-10005-70
pluriStrainer Mini, 10µm, 25pcs sterile in bag	43-10010-40
pluriStrainer Mini, 10µm, 50pcs, 25pcs sterile in bag	43-10010-50
pluriStrainer Mini, 10µm, 100pcs, 25pcs sterile in bag	43-10010-60
pluriStrainer Mini, 10µm, 500pcs, 25pcs sterile in bag	43-10010-70
pluriStrainer Mini, 20µm, 25pcs sterile in bag	43-10020-40
pluriStrainer Mini, 20µm, 50pcs, 25pcs sterile in bag	43-10020-50
pluriStrainer Mini, 20µm, 100pcs, 25pcs sterile in bag	43-10020-60
pluriStrainer Mini, 20µm, 500pcs, 25pcs sterile in bag	43-10020-70
pluriStrainer Mini, 40µm, 25pcs sterile in bag	43-10040-40
pluriStrainer Mini, 40µm, 50pcs, 25pcs sterile in bag	43-10040-50
pluriStrainer Mini, 40µm, 100pcs, 25pcs sterile in bag	43-10040-60
pluriStrainer Mini, 40µm, 500pcs, 25pcs sterile in bag	43-10040-70
pluriStrainer Mini, 70µm, 25pcs sterile in bag	43-10070-40
pluriStrainer Mini, 70µm, 50pcs, 25pcs sterile in bag	43-10070-50
pluriStrainer Mini, 70µm, 100pcs, 25pcs sterile in bag	43-10070-60
pluriStrainer Mini, 70µm, 500pcs, 25pcs sterile in bag	43-10070-70
pluriStrainer Mini, 100µm, 25pcs sterile in bag	43-10100-40
pluriStrainer Mini, 100µm, 50pcs, 25pcs sterile in bag	43-10100-50
pluriStrainer Mini, 100µm, 100pcs, 25pcs sterile in bag	43-10100-60
pluriStrainer Mini, 100µm, 500pcs, 25pcs sterile in bag	43-10100-70
pluriStrainer Mini, 200µm, 25pcs sterile in bag	43-10200-40
pluriStrainer Mini, 200µm, 50pcs, 25pcs sterile in bag	43-10200-50
pluriStrainer Mini, 200µm, 100pcs, 25pcs sterile in bag	43-10200-60
pluriStrainer Mini, 200µm, 500pcs, 25pcs sterile in bag	43-10200-70



過篩簡易多樣！孔徑選擇最多！